

“

**OMNILUX™, developed to harness the healing power of light in order to deliver effective treatments**

#### What is acne?

If you have acne, you're far from alone. By the age of 24, more than 85 percent of all adolescents and young adults will have had acne at one time or another – and up to 50% of all adults will experience acne after they turn 25.

Acne is one of the world's most common skin ailments and also one of the most misunderstood. It is becoming clear that diet or food items may not directly cause acne but, instead, acne is caused by an excess production of oil in the sebaceous glands – a natural response to the hormone imbalances common amongst young adults and pregnant women. This oil combines with skin cells to clog the glands. Normal skin bacteria colonise the pilo-sebaceous glands and the body's response causes the inflammation of reddening of the skin that we associate with acne.



Before



After (post 8 treatments)

There are numerous ways to treat acne. While acne cannot be prevented nor cured, researchers have identified a painless and effective treatment in the fight against acne.

**Talk to your skincare specialist and learn how OMNILUX™ can benefit you.**

“

*Realise the potential of light therapy*

Effective, Suitable for most skin types, Non-Invasive therapy, **NO side effects**, **NO pain**, **NO downtime**, Quick and simple, **Promotes body's own natural healing response**, Treatment régimes tailored specifically for you



**OMNILUX**  
leaders in light therapy

**ACNE**  
Mild to Moderate

#### Distributed by

True Solutions International, 15-17 Merton Street,  
Zetland NSW 2017, Australia  
Customer Service Australia 1800 808 993  
Customer Service New Zealand 0800 575 565

[www.spauniverse.com.au](http://www.spauniverse.com.au) [www.spauniverse.co.nz](http://www.spauniverse.co.nz)

The *healing* power of  
**LIGHT** therapy

“

*OMNILUX™ Light Therapy stimulates your body's natural processes to **heal** and **rejuvenate** the skin.*

#### What is OMNILUX™ Combination Light Therapy?


OMNILUX™ is a medical device that utilises LED technology. The OMNILUX™ Combination Light Therapy uses blue LEDs in conjunction with red LEDs to (1) destroy the bacteria responsible for the inflammation (2) stimulate cell growth (3) decrease inflammation and pore size and (4) stabilise oil production in the future.

#### How does it work?

One of the key bacteria responsible for the inflammation is Propionibacterium acnes or Pacnes. This bacteria produces natural chemicals called porphyrins, which are sensitive to light at different wavelengths. OMNILUX blue™ neutralises the bacteria that cause the redness and inflammation of acne. Combined with the anti-inflammatory properties of OMNILUX revive™, it helps minimise the redness of acne lesions and promotes healthier skin.

#### Who can benefit from the OMNILUX™ Combination Light Therapy treatment?

OMNILUX™ benefits anyone suffering with mild to moderate inflammatory acne. Consult your skincare specialist for a diagnosis of your acne. He or she will be able to tell you if you're a candidate.



“

*...benefits anyone suffering with mild to moderate inflammatory acne.*



#### What happens during treatment?

You will be instructed to remove any makeup or sunscreen before the treatment. Safety eyewear will be provided to you during treatment sessions as the light will be positioned close to your face. Each treatment session will take approximately 20 minutes. You will experience a feeling of relaxation and well-being during and after the treatment. Following your 20-minute session, you can go right back to your normal routine.

#### What results can be expected from a series of treatments?

Optimum results will be seen between 4 to 8 weeks after treatment. The light stimulates natural processes that go on after the treatment has stopped. Studies have shown complete clearance of lesions. On average you should expect 70-80% of acne lesions to clear.

#### How long do the effects last?

For most patients, acne is gone for 6 months. In some cases, touch-up treatments (typically 2 to 4 treatments) may be necessary.

“

*It's completely pain-free.*

**CLINICALLY PROVEN  
FOR PROBLEM SKIN.  
81% CLEARANCE  
REPORTED IN TRIALS.\***

#### How many treatments will I need?

Over the course of 4 weeks, 8 OMNILUX treatments will be administered, alternating between blue and red LEDs.

Depending on your type of acne, your skin specialist may also suggest using a face wash and/or a gel containing salicylic acid for optimum results.

Your skincare specialist can adjust the number of treatments if he or she find it appropriate. Consult your skincare specialist to determine your appropriate treatment program.

\*Goldberg DG, Russell B; Combination blue (415nm) and red (633 nm) LED phototherapy in the treatment of mild to severe acne vulgaris. Journal of Cosmetic and Laser Therapy. 2006; 8: 71-75.